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## 3M™ Sterilization Assurance

Technical Information

### Performance of the 3M™ Attest™ Rapid Readout Biological Indicator 1292 in 121°C Vacuum Assisted Steam Sterilization Cycles

#### Background

The Instructions for Use (IFU) provided with each box of 3M™ Attest™ Rapid Readout Biological Indicators 1292 states the product's FDA-cleared indications for use as follows:

- Use the Attest 1292 Rapid Readout Biological Indicator to monitor:
1. 250°F (121°C) gravity steam sterilization cycles.
  2. 270°F (132°C) vacuum assisted steam sterilization cycles.

As described in AAMI ST79, when monitoring steam sterilization cycles, biological indicators should be used within a representative Process Challenge Device (PCD). It is up to the user to determine the appropriate PCD for a particular sterilizer, such as for sterilizers >2 ft<sup>3</sup>, a biological indicator in a biological indicator challenge test pack (e.g., 3M™ Attest™ Rapid Readout Biological Indicators 1292 in a 16-towel PCD) or in a commercially available disposable PCD (e.g., the 3M™ Attest™ Steam Biological Test Pack 1296).

In response to customer inquiries about using this biological indicator to monitor 121°C vacuum assisted steam sterilization cycles, a study was conducted to validate the product's performance at this temperature. The study investigated key biological indicator performance requirements including D-value, Survival/Kill response, and Reduced Incubation Time.

#### Method

Three lots of 3M™ Attest™ Rapid Readout Biological Indicators 1292 were used. Sterilization cycles were performed in a steam resistometer (H&W Technologies, Rochester, NY) at 121°C. Exposed biological indicators were activated and incubated for 3 hours in a 3M™ Attest™ Auto-reader 390 to determine the 3-hour fluorescent result. The biological indicators were then transferred to a bacterial incubator set at 60°C for the pH color change incubation period.

D-value at 121°C was determined for each lot using the Holcomb-Spearman-Kärber method.

Survival/Kill times were calculated using the formulas provided in ISO 11138-1:2006. Fifty (50) biological indicators from each lot were exposed at 121°C for the calculated survival time in a steam resistometer. Another fifty (50) biological indicators from each lot were exposed at 121°C for the calculated kill time in a steam resistometer. The 3-hour fluorescent and 7-day visual pH color change responses were recorded for each biological indicator (BI).

Reduced Incubation Time was evaluated using the following method. Biological indicators from each lot were exposed to partial sterilization cycles at 121°C in a steam resistometer that resulted in between 30% and 80% growth positive results after 7 days of incubation. After exposure, the

biological indicators were activated and incubated in 3M™ Attest™ Auto-readers 390 for 3 hours and the fluorescence result was recorded. The biological indicators were then transferred to bacterial incubators and the visual pH color change response was recorded at 7 days.

#### Results

A summary of the validation acceptance criteria and results is provided in the table below. The acceptance criteria were satisfied for each performance attribute investigated.

		3M™ Attest™ BI 1292 Lot #			
		Acceptance Criteria	2016-07PB	2016-10TN	2016-12PB
Survival at Calculated Survival Time	D-value at 121°C	≥1.5 minutes	Accept	Accept	Accept
	3-hour Fluorescent Result	All BIs Positive	Accept	Accept	Accept
	7-day pH Color Change Result	All BIs Positive	Accept	Accept	Accept
Kill at Calculated Kill Time	3-hour Fluorescent Result	All BIs Negative	Accept	Accept	Accept
	7-day pH Color Change Result	All BIs Negative	Accept	Accept	Accept
Reduced Incubation Time	3-hour Fluorescent Result	≥97% Correlation to 7-day Growth Result	Accept	Accept	Accept

#### D-value

The D-value for each lot met the minimum 1.5 minutes as required in the FDA BI guidance document.

3M™ Attest™ BI 1292 Lot #	D-value (minutes)
2016-07PB	1.5
2016-10TN	1.8
2016-12PB	1.6

## Performance of the 3M™ Attest™ Rapid Readout Biological Indicator 1292 in 121°C Vacuum Assisted Steam Sterilization Cycles

### Survival/Kill

All BIs survived the calculated survival exposure times at 121°C and all BIs were negative after exposure to the calculated kill times at 121°C.

3M™ Attest™ BI 1292 Lot #			2016-07PB	2016-10TN	2016-12PB
Survival	Survival at Calculated Survival Time	# Positive 3-hour Fluorescent Result	50/50	50/50	50/50
		# Positive 7-day pH Color Change Result	50/50	50/50	50/50
		Calculated Survival Time min:sec	5:29	6:14	5:25
Kill	Kill at Calculated Kill Time	# Positive 3-hour Fluorescent Result	0/50	0/50	0/50
		# Positive 7-day pH Color Change Result	0/50	0/50	0/50
		Calculated Kill Time min:sec	14:29	17:02	15:01

### Reduced Incubation Time

Biological indicators from each lot were exposed to fractional 121°C sterilization cycles. The 3-hour results were compared to the 7-day growth results. All three lots satisfied the FDA's 97% acceptance criteria for reduced incubation time.

3M™ Attest™ BI 1292 Lot #			2016-07PB	2016-10TN	2016-12PB
# of Growth Positives at 7 days			108	217	84
Validation of Reduced Incubation Time		# of False Negatives at 3 Hours	0	2	0
		Sensitivity* at 3 hours (%)	100	99	100

\*Sensitivity =  $\frac{(\text{Number of Growth Positives after 7 days}) - (\text{Number of False Negatives})}{\text{Number of Growth Positives after 7 days}} \times 100$

### Conclusion

The acceptance criteria were satisfied for each performance attribute investigated in this study. It is therefore concluded that the 3M™ Attest™ Rapid Readout Biological Indicator 1292 is validated for use in monitoring 121°C vacuum assisted steam sterilization cycles.

### References

1. Data on file at 3M (TECH-REPORT-05-281984)
2. 3M™ Attest™ Rapid Readout Biological Indicators 1292, written IFU, © 2015
3. U.S. Food and Drug Administration: *Guidance for Industry and FDA Staff, Biological Indicator (BI) Premarket Notification [510(k)] Submissions*. October 4, 2007.
4. ISO: *Sterilization of health care products - Biological indicators*, ISO 11138-1,3:2006.
5. Association for the Advancement of Medical Instrumentation. *Comprehensive guide to steam sterilization and sterility assurance in health care facilities*. ANSI/AAMI ST79: 2010 and A1:2010 and A2:2011 and A3:2012 and A4:2013 (Consolidated Text).

**Infection Prevention Division**  
**3M Health Care**  
3M Center, 2510 Conway Ave  
St. Paul, MN 55144-1000  
1-800-228-3957  
www.3M.com/infectionprevention

Please recycle. Printed in USA.  
© 2016. All rights reserved.  
Issued: 3/16 11460HB  
70-9700-0044-7

3M and Attest are trademarks of 3M Company.  
Used under license by 3M subsidiaries and affiliates.